



**Triosite®
Bioactive
Ceramic Bone
Graft Substitute**



A Cost Effective Bioactive Approach to Bone Augmentation



zimmer | trauma
Confidence in your hands®



A Cost Effective Bioactive Approach to Bone Augmentation

In orthopaedics and traumatology, surgeons are often faced with reconstruction of lost, diseased or damaged bone. Traditionally, the use of autologous grafts has been considered ideal, however the supply of material is restricted, and harvesting autograft bone requires additional trauma and pain to the patient and associated blood loss. Allograft material may be available via local or regional bone banks, however its performance is considered to be inferior to autografts and there is potential for the material to act as a vector for microbiological contamination.

Synthetic Solution

There is a need for a biocompatible synthetic material that can be used to replace or augment graft material. Ideally, the material should become integrated with the surrounding bone tissue and ultimately be replaced by new healthy bone, as is seen with autologous graft.

Triosite Bone Graft Substitute is a bioactive calcium phosphate ceramic composed of hydroxyapatite (HA) and tricalcium phosphate (TCP). It has micro- and macropores and is supplied in both block and granular form. The material is soluble and gradually dissolves in the body, seeding new bone formation as it releases calcium and phosphate ions into the biological milieu. With time, the porous structure becomes completely infiltrated with, and replaced by healthy viable bone.

Triosite Bone Graft Substitute possesses the unique properties of controlled bioactivity and solubility as a result of its chemical composition. It consists of an optimum balance of the more stable phase of HA and more soluble TCP.

Triosite Bone Graft Substitute is indicated to augment or substitute for bone graft in nonloadbearing clinical applications.

Proven Clinical Benefits*

The effectiveness of *Triosite* Bone Graft Substitute as has been critically evaluated in a large prospective randomized study involving 341 patients undergoing posterior spinal fusion with associated instrumentation⁽¹⁾. The performance of *Triosite* Bone Graft Substitute was characterized by:

- Equivalent clinical outcome to autologous bone graft
- Fewer spinal wound healing problems than autologous grafts
- Avoidance of donor site pain and infections
- No allergic reactions

Histological analysis of biopsies containing *Triosite* Bone Graft Substitute show that it provides a favorable scaffolding for the stimulation and growth of new bone, and gradual incorporation into the fusion mass. Progressive *Triosite* Bone Graft Substitute resorption involved the same amount of well-differentiated bone ingrowth. It should be noted, however, that *Triosite* Bone Graft Substitute must be surrounded by and covered with host bone. As with other prosthetic materials, extensive contact with soft tissues can lead to fibrous encapsulation and delay or prevent full incorporation and replacement by bone.

* See publications at the end of this brochure.




Figure 1
Macropores (300–600 microns), which are a network of interconnected spaces and allow the biological infiltration and cellular colonization by osteoblasts and osteoclasts. (Data on file at Biomatlante)




Figure 2
Micropores (<10 microns) which are the intercrystalline spaces where dissolution and recrystallization take place. (Data on file at Biomatlante)

***Triosite* Bone Graft Substitute has also been used successfully in other clinical situations, such as the filling of pathogenic or mastoid cavity bone defects^(2, 3), and to replace bone in the treatment of long-bone fractures⁽⁴⁾.**

An Overview of Triosite Biphasic Ceramic Bone Substitute

Hydroxyapatite is stoichiometrically equivalent to the crystals of biological apatite found in natural bone. However, many chemical substitutions are found in biological apatite, such as calcium carbonates, which lead to imperfections in the crystal structure affecting size shape and particularly dissolution properties. Of all the synthetic calcium phosphate ceramics, hydroxyapatite is one of the least soluble and bioactive⁽⁵⁾.

In order to enhance the bioactivity of hydroxyapatite, it can be combined with tricalcium phosphate, which is soluble and reactive in biological terms. By controlling the relative proportions of these constituents it is possible to produce a ceramic with precise biological properties.

Based upon several studies^(6, 7, 8), Zimmer selected a balance of 60% HA and 40% TCP for *Triosite* Bone Graft Substitute because it appears to provide the optimum performance in terms of controlled resorption and osseous substitution.

The porosity of the material also has a marked influence on its solubility and bioactivity⁽⁹⁾, so the physical properties of *Triosite* Bone Graft Substitute are also carefully controlled. The pores in *Triosite* Bone Graft Substitute can be divided into two categories: the micropores (Fig. 2) and the macropores (Fig. 1).

At the microscopic level, biological fluids seep through micropores, and become enriched with calcium and phosphate ions which are released during dissolution of TCP. As the solution becomes saturated, crystals of biological apatite precipitate in the spaces and bond the substrate together^(10, 11, 12) integrating non-collagenic proteins like growth factors and osteogenic molecules. An increase in the mechanical properties of the biomaterial/tissue construct is apparent⁽¹³⁾.

At the histological level, viable lamellar bone (O) is deposited directly onto the surfaces of the ceramic scaffold (Fig. 3). Over time the lamellar bone is remodelled to Haversian bone (H), which also gradually replaces the *Triosite* Bone Graft Substitute particles (T) (Fig. 4).^(14, 15)




Figure 3
Lamella bone (O, area in pink)




Figure 4
Haversian bone (H)
Triosite Bone Graft Substitute particles (T)

Technical Properties

Triosite Bone Graft Substitute is a high-purity ceramic which is synthesized using a special process involving high-temperature sintering.

Typical Composition	60% Hydroxyapatite ($\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2$) 40% Tricalcium Phosphate ($\text{Ca}_3(\text{PO}_4)_2$)
Ca:P ratio (by mass)	1.6
Purity	< 50 ppm Heavy Metals
Porosity	60–70%
Pore Size	Macropores = 300–600 microns Micropores \leq 10 microns

References

1. Ransford, A.O., et al., "Synthetic porous ceramic compared with autograft in scoliosis surgery. A prospective, randomized study of 341 patients." *J Bone Joint Surg Br*, 1998. 80(1): 13-8.
2. Gouin, F., Delecrin, J., Passuti, N., Touchais, S., Poirier, P., Bainvel, J. "Biphasic macroporous calcium phosphate ceramic bone substitute for filling bone defects: A report of 23 cases." *Revue de Chirurgie Orthopedique*: 81; 59–65 (1995).
3. Daculsi, G., Corlieu, P., D'Arc, M., Gersdorff, M. "Macroporous biphasic calcium phosphate efficiency in mastoid cavity obliteration: experimental and clinical findings." *Annals of Otolaryngology, Rhinology & Laryngology*: 101 (8); 669–674 (1992).
4. Daculsi, G., Passuti, N., Martin, S., Deudon, C., Legeros, R., Raheer, S. "Macroporous calcium phosphate ceramic for long bone surgery in humans and dogs. Clinical and histological study." *Journal of Biomedical Materials Research*: 24; 379–396 (1990).
5. Legeros, R., Parsons, J., Daculsi, G., Driessens, F., Lee, D., Liu, S., Megsger, D., Peterson, D., Walker, M. "Significance of the porosity and physicochemical chemistry of calcium phosphate ceramics. Biodegradation and bioresorption." *Ann. N.Y. Acad. Sci.*: 523; 268–271 (1998).
6. Daculsi, G., Legeros, R., Nery, E., Lynch, K., Kerbel, B. "Transformation of biphasic calcium phosphate ceramics in vivo: Ultrastructural and physicochemical characterisation." *J. Biomed. Mater. Res.*: 23; 883–894 (1989).
7. Nery, E., Lynch, K., Hirthe, W., Mueller, K. "Bioceramics implants in surgically produced infrabony defects." *J. Periodontol.*: 46; 328 (1975).
8. Nery, E., Legeros, R., Lynch, E., Kalbleish, . "Tissue response to biphasic ceramic with different ratios of sintered apatites." *J. Dent. Res.*: 65; 292(1986).
9. Daculsi, G., Passuti, N. "Effect of macroporosity for osseous substitution of calcium phosphate ceramics." *Biomaterials*: 11; 86–87 (1990)
10. Daculsi, G., Legeros, R., Heughebaert, M., Barbieux, I. "Formation of carbonate apatite crystals after implantation of calcium phosphate ceramics." *Calcif. Tissue Int.*: 46; 20–27 (1990).
11. Daculsi, G., Legeros, R., Mitre, D. "Crystal dissolution of biological and ceramic apatites". *Calcif. Tissue Int.*: 45; 95–103 (1989).
12. Daculsi, G., Legeros, R., Deudon, C. "Scanning and transmission electron microscopy and electron probe analysis of the interface between implants and host bone. Osseo-coalescence versus osseo-integration." *Scann. Microscopy*: 4; 309–314 (1990).
13. Trecant, M., Delecrin, J., Nguyen, J., Royer, J., Passuti, N., Daculsi, G. "Influence of post-implantation physico-chemical changes in a macroporous ceramic on its mechanical strength." *J. Mater. Science: materials in medicine*: 7; 227–229 (1996).
14. Passuti, N., Daculsi, G., Rogez, J., Martin, S., Bainvel, J. "Macroporous calcium phosphate ceramic performance in human spine." *Clin. Orthop. Relat. Res.*: 248; 169–176 (1989).
15. Passuti, N., Daculsi, G., Martin, S., Deudon, C. "Macroporous polycrystalline calcium phosphate implant for spinal fusion in man and dogs." *In Handbook of Bioactive Ceramics, vol. 2, CRC Press*; 345–354 (1991).

Disclaimer

This document is intended exclusively for physicians and is not intended for laypersons.

Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part.

Information contained in this document was gathered and compiled by medical experts and qualified Zimmer personnel. The information contained herein is accurate to the best knowledge of Zimmer and of those experts and personnel involved in its compilation. However, Zimmer does not assume any liability for the accuracy, completeness or quality of the information in this document, and Zimmer is not liable for any losses, tangible or intangible, that may be caused by the use of this information.



Triosite Sticks

Size (mm)	Quantity	REF
5.0 × 5.0 × 10.0	Box of 2	00-1109-113-01
5.0 × 5.0 × 10.0	Box of 4	00-1109-113-02
5.0 × 5.0 × 20.0	Box of 2	00-1109-114-01
5.0 × 5.0 × 20.0	Box of 4	00-1109-114-02



Triosite Cones

Size (mm)	Quantity	REF
8.0	1	00-1109-085-15
12.0	1	00-1109-132-15



Triosite Osteotomy Wedges

Size (mm)	Quantity	REF
3.0 × 6.0*	1	00-1109-000-36
4.0 × 8.0*	1	00-1109-000-48
5.0 × 10.0*	1	00-1109-005-10
6.0 × 12.0*	1	00-1109-006-12



Triosite Granules

Size (mm)	Quantity	REF
0.5–1.0	2 cc	00-0060-130-00
2.0–3.0	10 cc	00-1109-115-01
2.0–4.0	16 cc	00-1109-115-02
2.0–3.0	5 cc	00-1109-005-03



Triosite Disks

Size (mm)	Quantity	REF
25.0 × 5.0	1	00-1109-005-25

* Note: Wedges are 30 mm × 30 mm square. Dimensions quoted in the product description relate to minimum and maximum thickness; the latter is a good approximation of the wedge angle.

CE 0123

Manufacturer: Biomatlante
Distributor: Zimmer GmbH

Contact your Zimmer representative or visit us at www.zimmer.com

