

# Pocket Reference Chart





## Standard Frames

Cat. No.	Description
7107-1442	Medium struts, 155mm rings
7107-1551	Short struts, 180mm rings
7107-1552	Medium struts, 180mm rings



## Open Frames

Cat. No.	Description
7107-2552	Medium struts, 180mm rings
7107-2662	Medium struts, 205mm rings




## Foot Frame

Cat. No.	Description
7107-3442	Medium struts, 155mm rings



## Full Rings

Cat. No.	Inner Diameter
7107-0111	80mm
7107-0112	105mm
7107-0113	130mm
7107-0114	155mm
7107-0115	180mm
7107-0116	205mm
7107-0117	230mm



## Half Rings

Cat. No.	Inner Diameter
7107-0123	130mm
7107-0124	155mm
7107-0125	180mm
7107-0126	205mm
7107-0127	230mm



## 2/3 Rings

Cat. No.	Inner Diameter
7107-0131	80mm
7107-0132	105mm
7107-0133	130mm
7107-0134	155mm
7107-0135	180mm
7107-0136	205mm
7107-0137	230mm



## Foot Rings & U-Plates

Cat. No.	Description
7107-0141	80mm
7107-0142	105mm
7107-0143	130mm
7107-0144	155mm Long
7107-0145	180mm Long
7107-0154	155mm Short
7107-0155	180mm Short



## Accessory Block

Cat. No.	Description
7107-0161	Accessory Block



## Standard Struts

Cat. No.	Description
7107-0200	XX-Short, 59-76mm
7107-0205	X-Short, 75-96mm
7107-0210	Short, 90-125mm
7107-0220	Medium, 116-178mm
7107-0230	Long, 169-283mm



## Lobe

Cat. No.	Description
7107-0160	



## FAST FX<sup>o</sup> Struts

Cat. No.	Description
7107-0705	X-Short, 91-121mm
7107-0710	Short, 116-152mm
7107-0720	Medium, 143-205mm
7107-0730	Long, 195-311mm



## Bolts

Cat. No.	Description
7107-0310	Shoulder Bolt
7107-0515	Connection Bolt, 12mm



## Wrench

Cat. No.	Description
7107-0410	Torque Wrench




## ID Bands

Cat. No.	Description
7107-0320	Standard Strut ID Band Kit
7107-0340	FAST F <sub>x</sub> Strut ID Band Kit




## Trays & Caddies


Cat. No.	Description
7107-0490	Rings Tray



7107-0491	Struts and Wrenches Tray
-----------	--------------------------



7107-0492	Pins and Wires Tray
-----------	---------------------



7107-0493	Lid Insert
-----------	------------



7107-0494	Nuts, Bolts and Washers Caddy
-----------	-------------------------------



## Basic Tibial Set 7107-0500

---

### Bolts and Nuts

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-0600	Wire Fixation Bolt, Cannulated	6
10-0615	Half Pin Fixation Bolt, 5mm	6
10-0616	Half Pin Fixation Bolt, 6mm	6
10-0700	Wire Fixation Bolt, Slotted	12
10-3300	Nut, 10mm	40
7107-0310	Shoulder Bolt	24
7107-0515	Connection Bolt, 12mm	25

---

### Sockets

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-0910	Threaded Socket, 30mm	8
10-0911	Threaded Socket, 60mm	8

---

### Wires

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-2102	Smooth Wire, 1.8mm Diameter, 370mm	12
10-2107	Olive Wire, 1.8mm Diameter, 400mm	12

---

### Threaded Rods

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-2300	60mm	4
10-2301	80mm	4
10-2302	100mm	4
10-2303	120mm	4
10-2304	150mm	4
10-2305	200mm	4

---

### Washers

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-2600	Slotted, 4mm X 14mm	20
10-2701	2mm X 14mm	20

---

### Drills

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-3044	2.7mm with 6mm Shank	1
10-3045	3.8mm with 6mm Shank	1
10-3046	4.8mm with 6mm Shank	1

---

### Centering Sleeves

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-3404	4mm	6
10-3405	5mm	6
10-3406	6mm	6

---

### Rancho Cubes

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-3451	1-Hole	4
10-3452	2-Hole	4
10-3453	3-Hole	4
10-3454	4-Hole	4
10-3455	5-Hole	4
11-2727	Set Screw	20

## Titanium Half Pins

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
12-2704	4.0mm X 40mm	6
12-2710	5.0mm X 35mm	6
12-2727	6.0mm X 50mm	6

---

## Instruments

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
10-2900	Open/Box end, 10mm Wrench	2
10-3000	Tubular Angulated, 10 mm Wrench	1
10-3101	Dynamometric Wire Tensioner	1
11-2716	Pin Driver/Extractor	1
11-2717	Hex Driver, 90°	1
11-2719	Hex Driver, Straight	1
21-0064	Super-Bite Pin Cutter	1
7101-1005	Rancho Adapter Sleeve	2
7103-1040	Drill Guide with Trocar	2
7107-0410	Torque Wrench	1
L20-2002	10mm X 13mm Wrench	1

---

## Full Rings

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
7107-0114	155mm	2
7107-0115	180mm	2

---

## Half Rings

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
7107-0124	155mm	4
7107-0125	180mm	4

---

## 2/3 Rings

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
7107-0134	155mm	1
7107-0135	180mm	1

---

## FAST FX<sup>o</sup> Struts

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
7107-0705	X-Short	6
7107-0710	Short	6
7107-0720	Medium	6
7107-0730	Long	6

---

## ID Bands

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
7107-0340	FAST F <sub>x</sub> Strut ID Band Kit	3

---

## Trays and Cases

<b>Cat. No.</b>	<b>Description</b>	<b>Qty.</b>
7107-0490	Rings Tray	1
7107-0491	Struts and Wrenches Tray	1
7107-0492	Pins and Wires Tray	1
7107-0493	Lid Insert	1
7107-0494	Nuts, Bolts and Washers Caddy	1
7112-9400	Large Case, 4.8"	1
7112-9401	Small Case, 2.4"	1
7112-9402	Lid	2



---

**Orthopaedics** - Smith & Nephew, Inc., 1450 Brooks Road,  
Memphis, TN 38116, USA  
Phone: 1-901-396-2121, Information: 1-800-821-5700,  
Orders/inquiries: 1-800-238-7538  
[www.smith-nephew.com/orthopaedics](http://www.smith-nephew.com/orthopaedics)

Pictures not to scale.

30036215001d 7108-0452 03/06

°Trademark of Smith & Nephew. Reg. US Pat. & TM Off.